## RiO5 METHOD (10)

IRSN/LMRE/ MESURE/BETA/MO4010.11.1 Indice : 9

Sabine Charmasson

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Béatrice Boulet ([beatrice.boulet@irsn.fr](mailto:beatrice.boulet@irsn.fr))

Carole Vivien ([carole.vivien@irsn.fr](mailto:carole.vivien@irsn.fr))

# **$^{90}\text{Sr}$ — Chemical treatment of a sample for measurement of $^{90}\text{Sr}$ by proportional counter —*sediment and biota samples*—**

### **Disclaimer**

It is the responsibility of each analyst to follow established practices when handling and examining the samples referenced in this Rio5 Cookbook. Although the methods may have been tested by each laboratory identified as the source, each user must perform a validation procedure to ensure the validity of their results. Woods Hole Oceanographic Institution, its officers, directors and employees are not responsible for any of the data or the results that may be achieved from using the information in the Rio5 Cookbook and disclaim all liability for the same.

## Table of Contents

<b>1</b>	<b>SCOPE</b>	<b>1</b>
<b>2</b>	<b>EQUIPMENT AND CHEMICAL REAGENTS</b>	<b>1</b>
2.1	EQUIPMENT	1
2.2	TRACERS	1
2.3	CHEMICAL REAGENTS	2
2.4	SOLUTIONS	2
<b>3</b>	<b>PROCEDURE</b>	<b>3</b>
3.1	ÉTAPES DE RADIOCHIMIE/ <i>RADIOCHEMISTRY STEPS</i>	3
3.2	DETAIL DES ETAPES DE PREPARATION DES ECHANTILLONS	4
<b>4</b>	<b>REFERENCES</b>	<b>10</b>
<b>5</b>	<b>FLOW CHART</b>	<b>11</b>

## 1 SCOPE

This method specifies the minimum requirements and laboratory methods for the analysis of  $^{90}\text{Sr}$  in sediment or biota samples. It gives a description of the radiochemical steps for quantification of strontium 90 by gas proportional counter.

$^{90}\text{Sr}$  is measured in the samples of the environment by  $\beta$  counting of its descendant  $^{90}\text{Y}$ , with a proportional counter. The procedure requires prior radiochemistry to extract  $^{90}\text{Sr}$  from the sample and separate it from  $^{90}\text{Y}$  initially present. The measurement is carried out after about 14 days, after which yttrium 90 is at secular equilibrium with strontium 90.

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$^{90}\text{Sr}$  (at equilibrium with  $^{90}\text{Y}$ ), dilution of SR90ELSB30 from Areva CERCA – LEA (~ 50 Bq.g<sup>-1</sup>)

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## Reference material

IAEA-156 (clover), IAEA Soil-6, IAEA-321 (milk powder), IAEA-375 (soil), IAEA-152 (milk powder)

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*The purity of the reagents is verified by the analysis of the white sample.*

### 3 PROCEDURE

#### Warning

The original documents quote various internal documents relating in particular to the calibration, the maintenance of the equipment used as well as aspects of internal quality assurance in the laboratory. These documents are not published here and it is appropriate for each laboratory that uses these procedures to verify that it works in the state of the art.

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*The different steps, schematized on figure F1, are :*

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*The solution then obtained is analyzed to evaluate the chemical yield by measuring strontium 85 using gamma spectrometry or, stable strontium by elementary measurement. After a fortnight (at least) strontium 90 is in radioactive equilibrium with its son yttrium 90.*

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*The yttrium hydroxide is extracted from the solution by precipitation. This precipitate is dissolved again and the yttrium oxalate is precipitated and deposited on a filter.*

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*Counting of the filter by gas proportional counter.*

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*The procedure is detailed below (3.2). The glassware is specified for information only and for 20 g of ash.*

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*Weigh approximately 20 g of previously calcinated sample (480 ° C) - and if necessary dried to constant weight - into a 250 or 400 ml beaker. Note the weight*

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*Wet the ashes with ultra-pure water (getting a cake).*

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##### *Sample spiking with $^{85}\text{Sr}$ radioactive tracer*

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*Add by weigh the desired amount of  $^{85}\text{Sr}$  tracer solution to the wet ashes*

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*This step consists in dissolving the minerals present in the form of oxides in the matrix*

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*Add about 30 mL of concentrated  $\text{HNO}_3$ .*

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*Heat at a gentle temperature (80 °C) all day long by adding 5 to 10 ml of  $\text{H}_2\text{O}_2$  regularly in fractions of 1 mL. Shake regularly. Add 9M  $\text{HCl}$  or concentrated  $\text{HCl}$  if necessary.*

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*Evaporate to dryness. Repeat once from step B.1.*

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*Dissolve the dry residue with 20 to 50 ml 4M  $\text{HCl}$ . Bring to a boil for at least 10 minutes.*

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*Remove from the hotplate. Add 20 to 50 mL of ultra-pure H<sub>2</sub>O and allow to cool.*

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*Transfer the filtrate to a 1 or 2 L beaker.*

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*Add Sr(NO<sub>3</sub>)<sub>2</sub> solution to obtain approximately 10 mg stable Sr.*

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*Add 10 mL of 25% NH<sub>4</sub>Cl, a tip of spatula of citric acid and an amount of oxalic acid equivalent to the previously dissolved hot sample.*

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*Adjust the pH of the sample to 4.5 with concentrated NH<sub>4</sub>OH using the pH meter. Towards pH = 3, there is appearance of a white precipitate.*

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*Stir the solution for 1 hour while heating at 70-80 °C.*

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*Leave to settle outside the plate having removed the magnetic stirring rod.*

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*Recover the precipitate, and dry in an oven (about 12 hours at 80-100 °C).*

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*On leaving the furnace, let the samples cool in a desiccator. Weigh them*

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***For samples known to be rich in phosphates (milk, water, fish, cheese), go directly to step F.***

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*Resume the sample (which is in carbonate form) with 8M HNO<sub>3</sub> or concentrated.*

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*Heat to remove H<sub>2</sub>O<sub>2</sub> until the end of the effervescence.*

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*Extend the volume of the solution to 350 mL with ultra-pure H<sub>2</sub>O.*

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*Adjust the pH to 8.5 with concentrated NH<sub>4</sub>OH using the pH meter.*

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*Set aside the supernatant in a 1 L beaker.*

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*In the supernatant add 10 mL of 25%  $\text{NH}_4\text{Cl}$ , a tip of spatula of citric acid and 10 g of oxalic acid previously dissolved while heating.*

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*Adjust the pH to 4.5 with concentrated  $\text{NH}_4\text{OH}$  using the pH meter. Towards pH = 3, a white precipitate appears.*

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*Collect the precipitate and dry it in an oven (about 12 hours at 80-100 °C)*

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*On leaving the oven, let the samples cool in a desiccator.*

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*Dissolve the precipitate with not more than 30 mL 0.5 M  $\text{Al}(\text{NO}_3)_3$ , 6 M  $\text{HNO}_3$ .*

**F.2** Assembler la colonne (robinet, filtre, corps, entonnoir).

*Assemble the column (tap, filter, body, funnel).*

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*Condition the column with 20 mL 6M  $\text{HNO}_3$ .*

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*Load the sample at the top of the column*

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*Elute the strontium with 40 mL of a 0.05 M HNO<sub>3</sub> solution in a vial for measurement by gamma spectrometry*

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*The samples must undergo a gamma spectrometry for the measurement of <sup>85</sup>Sr. Wait at least 15 days for strontium 90 and yttrium 90 to be at radioactive equilibrium before proceeding with radiochemistry in step G.*

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*After the equilibrium standby, add approximately 10 mg of yttrium to a 100 mL glass centrifuge tube by weighing.*

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*Pour the solution from the gamma counting vial into the centrifuge tube.*

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*Add NH<sub>4</sub>OH until a white precipitate of yttrium hydroxide appears (check with pH paper that the pH is between 8 and 9). **Note the date and time***

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*Resume the precipitate with a few drops of 6 M or 8 M  $\text{HNO}_3$ .*

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*Extend the solution to a volume between 30 and 50 mL with ultra-pure  $\text{H}_2\text{O}$ .*

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*Add 6 mL of a saturated solution of  $\text{H}_2\text{C}_2\text{O}_4$  until the appearance of a white yttrium oxalate co-precipitate.*

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*Heat for 20 minutes in a water bath set at 80°C.*

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*Place a pre-weighed cellulose acetate filter on the filter-holder of the filtration system. Wet the filter with milliQ water. Using a plastic pipette, take the coprecipitate and carefully place it in the center of the filter.*

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*Rinse the filter with milliQ water and ethanol.*

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*Carefully remove the filter put it in an oven (105 °C maximum) for 30 minutes minimum.*

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*Start counting immediately. **Note the date and time.***

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Standard NF ISO 18589-5 (may 2009) Measurements of radioactivity in the environment - Soil - Part 5 : measurement of strontium 90 -

IAEA N° 295, Measurement of radionuclides in food and the environment. A guidebook. Technical reports series, 1989.

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
Figure F1. Schematic representation of the steps for the preparation of yttrium oxalate