## RiO5 METHOD (11)

IRSN-LMRE/MESURE/ALPHA/MO4010.9.1

Sabine Charmasson

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#### **Disclaimer**

It is the responsibility of each analyst to follow established practices when handling and examining the samples referenced in this Rio5 Cookbook. Although the methods may have been tested by each laboratory identified as the source, each user must perform a validation procedure to ensure the validity of their results. Woods Hole Oceanographic Institution, its officers, directors and employees are not responsible for any of the data or the results that may be achieved from using the information in the Rio5 Cookbook and disclaim all liability for the same.

## Table of Contents

<b>1</b>	<b>SCOPE</b>	<b>1</b>
<b>2</b>	<b>EQUIPMENT CHEMICAL REAGENTS</b>	<b>1</b>
<b>2.1</b>	<b>EQUIPMENT</b>	<b>1</b>
<b>2.2</b>	<b>TRACERS</b>	<b>2</b>
<b>2.3</b>	<b>CHEMICAL REAGENTS</b>	<b>2</b>
<b>2.4</b>	<b>SOLUTIONS</b>	<b>2</b>
<b>3</b>	<b>PROCEDURE</b>	<b>3</b>
	<b>ANNEX 1. PREPARATION DES COLONNES/ <i>COLUMN PREPARATION</i></b>	<b>27</b>
	<b>ANNEX 2. PRESENCE D'INTERFERENCES SPECTRALES</b>	<b>29</b>
<b>4</b>	<b>REFERENCES</b>	<b>29</b>
<b>5</b>	<b>FLOW CHART</b>	<b>30</b>

## 1 SCOPE

This method both in French (original language) and *English (in italic characters)* specifies the minimum requirements and laboratory methods for the analysis of  $^{238}\text{Pu}$ ,  $^{239+240}\text{Pu}$ ,  $^{241}\text{Am}$  and  $^{244}\text{Cm}$  in sediment/soil and biota samples.

This procedure has two variants:

- procedure for iron-rich samples, for example : sediment (test portion (tp) between 100 g and 200 g or according to expected activity), soils (tp 100 g to 200 g), grassland, bryophytes, algae (tp 50 g)

- procedure for samples rich in phosphate ions, for example: fish, molluscs, crustaceans, dairy products, foodstuffs (tp 50 g)

After collection in the field, the samples are either frozen or keep fresh until return to the laboratory where they are dried at 80°C or freeze-dried) (they can be dried or freeze-dried directly on board) then ashed at 480°C (by increasing the furnace temperature step by step). Then the dissolution of the inorganic oxides is carried out followed by three co-precipitations allowing the actinides to be concentrated. A first load of the solution on column makes it possible to separate the fraction containing Pu from that containing Am and Cm. The two fractions are then treated according to different protocols. Various column loads make it possible to purify, on the one hand, the Pu fraction and, on the other hand, the Am-Cm fraction. Then, irrespective of the radionuclide to be measured, the purified fraction is electroplated and counted by alpha spectrometry.

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See text of Procedure for other materials

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*The use of the different devices is explained in the associated manufacturer's procedures.*

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*The laboratory glassware will be annotated ad minima with the sample number.*

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*See text of Procedure for other chemical reagents (e.g. resins).*

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*The quality of the reagents is checked by the analysis of the "BLANK".*

## 2.4 Solutions

- Concentrated nitric acid or 8 M  $\text{HNO}_3$
- 27.5%  $\text{FeCl}_3$  iron chloride
- 1M  $\text{HNO}_3$  / 93%  $\text{CH}_3\text{OH}$
- 0.1M  $\text{HCl}$  / 0.5M  $\text{NH}_4\text{SCN}$  / 80%  $\text{CH}_3\text{OH}$
- Concentrated  $\text{HCl}$ , 9 M  $\text{HCl}$ , 12 M  $\text{HCl}$
- Aluminum nitrate 1M  $\text{Al}(\text{NO}_3)_3$
- 0.1 M  $\text{MnCl}_2$
- 0.2 M  $\text{KMnO}_4$ .

See text of Procedure for other solutions.

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*after having weighed the ashed sample (step A), the dissolution of the inorganic oxides (step B) is followed by three co-precipitations (steps C / C', D / D', E) allowing the actinides to be concentrated (steps C and D for the iron-rich samples, steps C' and D' for those rich in phosphate ions,*

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*a first load of the solution on column makes it possible to separate the fraction containing Pu from that containing Am and Cm (step F),*

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*the two fractions are then treated according to different protocols. Various column loads make it possible to purify, on the one hand, the Pu fraction (step G) and, on the other hand, the Am-Cm fraction (steps H, I, J, K)*

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*then, irrespective of the radionuclide to be measured, the preparation of the source by electroplating (step L) and counting (step M) is carried out.*

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*The precisions on the glassware and the volumes are given for information only and for a maximum test portion (200 g for the soils and sediments, 50 g for the other matrices). The choice of glassware (size and shape) is left to the discretion of the operator.*

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*For small test portions of iron-rich samples (<10 g), if in stage C the quantity of calcium oxalate is not too great, one can proceed directly to step F*

**Warning**

The original documents quote various internal documents relating in particular to the calibration, the maintenance of the equipment used as well as aspects of internal quality assurance in the laboratory. These documents are not published here and it is appropriate for each laboratory that uses these procedures to verify that it works in the state of the art.

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***Addition of radioactive tracers <sup>242</sup>Pu and <sup>243</sup>Am (spiking)***

- A.7** Mettre des gants. *Put on gloves*
- A.8** Disposer le bécher contenant les cendres et le flacon contenant la solution de traçage sous la hotte.  
*Place the beaker containing the ashes and the vial containing the tracing solution under the hood.*
- A.9** Ouvrir le flacon de la solution de traçage avec précaution.  
*Open the vial of the tracer solution carefully.*

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***Addition of tracer volume greater than 0.5 mL (spiking)***

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**Add by weighing a known quantity of tracer (<0.5ml) (spiking)**

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*Print the sample number. Pour on the ashes the desired quantity by weight of tracer, knowing that a drop weighs approximately 40 mg (33 µL).*

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*Repeat these last two operations for each sample in the group.*

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*Pour the excess solution into the vial. Discard the dropper bottle into the radioactive waste bin.*

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*Paste the label printed by the scale (or note the weigh) onto the weighing tracking sheet.*

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*This step consists of dissolving the minerals present as oxides in the matrix. This operation does not make it possible to dissolve the silica which constitutes the major part of the insoluble fraction (steps B1 to B16). However, a small fraction of silica can be entrained during the washing of the ashes. This embarrassing silica for the rest of the procedure must be eliminated (steps B17 to B24).*

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*For some specific studies (plants, foodstuffs) it may be necessary to dissolve the entire insoluble fraction with hydrofluoric acid (steps B25 to B32)*

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*The volumes of reagents indicated correspond to the average volumes. The amount of reagent is also related to the amount of the initial test fraction.*

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*Add the following concentrated nitric acid volumes:*

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*Gently introduce the first 50 to 100 mL in 10 mL fractions to avoid splashes on the walls. Between each addition of HNO<sub>3</sub>, cover with a watch glass.*

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*Add 50 mL (for sols, sediments) / 20 mL (for other matrices) of H<sub>2</sub>O<sub>2</sub> per 10 mL fraction and leave to act until the end of the effervescence. Between each addition of H<sub>2</sub>O<sub>2</sub>, cover with a watch glass.*

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*Heat at a mild temperature (80 °C) for at least 8 hours, adding H<sub>2</sub>O<sub>2</sub> evenly in 10 to 30 mL fractions. Shake regularly. If the rust color of the sample persists add 50 to 100 mL of concentrated HCl.*

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*Transfer the supernatant into a 1 L [600 mL] beaker, reduce the volume to crystallization limit.*

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*Dissolve the pellet with 8M or concentrate HNO<sub>3</sub>:*

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*Heat gently (at 80 °C) for 8 hours. Agitate and add regularly a few tens of mL of H<sub>2</sub>O<sub>2</sub>.*

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*Add the supernatant obtained to the previous supernatant stored in the 1 L beaker*

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*Wash the insoluble fraction with 200 mL of ultra-pure water.*

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*Centrifuge and join the supernatant to the others already preserved.*

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*The insoluble fraction is discarded.*

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*Reduce the volume of the supernatant until flocculation of the silica, ie down to 100-300 mL*

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*Reserve the supernatant in a beaker.*

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*Wash the centrifugation pellet with 50 to 100 mL of ultra-pure water.*

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*In a Teflon beaker, transfer the centrifugation pellet by resuspending the silica in 5 to 10 mL of 8 M or concentrated HNO<sub>3</sub>.*

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*Bring to dryness on a sand bath (insert a stainless steel plate so that the beaker does not melt).*

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*Add a few tens of mL of HF and bring to dryness.*

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*Dissolve the residue with 20 mL of concentrated HNO<sub>3</sub> and 5 to 20 mL of H<sub>2</sub>O<sub>2</sub>. If residue remains, filter through cellulose acetate filter and discard the residue.*

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*Pour the recovery solution or the filtrate into the beaker.*

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***Dissolution of the « insoluble fraction » in the case of particular applications for certain studies (plants, foodstuffs).***

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***Precipitation of calcium oxalate (Iron-rich sample)***

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*Adjust the pH to 1.5 with concentrated  $\text{NH}_4\text{OH}$ . If the pH of the solution is initially less than -0.3, add concentrated ammonia without the titrator to increase the pH to about 0.5. Towards  $\text{pH} = 1$ , a white precipitate appears. If the precipitate is not visible, add 5 to 10 mL of the  $\text{CaCl}_2$  solution ( $100 \text{ mg.mL}^{-1}$ ).*

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*Dry the precipitate in an oven for at least 24 hours. After drying, transfer the precipitate into a beaker of sufficient capacity (minimum 1 L) previously tared*

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*Cover with a watch glass. Calcinate at 550°C (with increase of temperature step by step)*

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*Right out of the oven, let the samples cool down in the form of carbonate. Weigh them. Post the value on the Weighing Tracking Sheet.*

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### ***Precipitation of ferric hydroxide (iron-rich sample)***

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*Wet the carbonate (C.22) with ultrapure water.*

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*Dissolve the precipitate in a fraction of 10 mL of concentrated or 9 M HCl until complete dissolution. An effervescence occurs. Cover with a watch glass.*

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*Add 10 to 20 mL of H<sub>2</sub>O<sub>2</sub> and heat until the end of the effervescence.*

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*If a black color persists, filter the sample on a cellulose acetate filter and calcine it. Dissolve the calcination residue with 20 mL of 9M HCl and transfer to the beaker.*

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*Adjust with ultrapure water the volume of the solution at  $\frac{3}{4}$  the height of the beaker*

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*Adjust to pH 8.5 with concentrated NH<sub>4</sub>OH using the previously calibrated titrator. A rusty brown colloidal precipitate is formed. If the precipitate is not visible, add 0.5 to*

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***Precipitation of manganese oxide (sample rich in phosphate ions)***

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*Add 5 mL of 0.1 M  $\text{MnCl}_2$  and 1-2 mL of 0.2 M  $\text{KMnO}_4$ . After the brown precipitate of  $\text{MnO}_2$  appears, the solution is brought back to acidic pH (about 1.5 - 2) until it is clear. If the brown precipitate does not appear before acidification, the solution may still contain  $\text{H}_2\text{O}_2$ : heat the solution and continue the precipitation without adding  $\text{MnCl}_2$ .*

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*The pH is manually adjusted with  $\text{NH}_4\text{OH}$  to pH 3.5 (or up to 4 if the amount of phosphate ions is low) up to the phosphate precipitation limit (white veil).*

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***3<sup>rd</sup> precipitation of manganese oxide (sample rich in phosphate ions)***

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*If too much phosphate has been co-precipitated, restart the precipitation from step C'.6 without adding  $\text{MnCl}_2$ .*

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*Dissolve the precipitates with  $\text{HNO}_3$  and a few mL of  $\text{H}_2\text{O}_2$ .*

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*Heat until the end of the effervescence of  $\text{H}_2\text{O}_2$  to eliminate it completely from the solution (otherwise the precipitate will not form).*

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*Transfer to a 2 L beaker and adjust the volume of the solution to 1500 mL with ultrapure water.*

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*Adjust the pH to 8 with concentrated  $\text{NH}_4\text{OH}$ . Close to pH = 3 - 4, a white precipitate of calcium phosphate appears.*

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*Add 2 to 3 mL of 0.2 M  $\text{KMnO}_4$ . After the brown precipitate of  $\text{MnO}_2$  appears, the solution is brought back to acidic pH (about 1.5 - 2) until it is clear. The pH is manually adjusted with  $\text{NH}_4\text{OH}$  to pH 3.5 (or up to 4 if the amount of phosphate ions is low) up to the phosphate precipitation limit (white veil).*

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*Leave to settle. Withdraw the supernatant.*

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*Combine the supernatants and throw them away.*

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*Dissolve 20g of oxalic acid in 200 ml of water while heating.*

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*Dissolve the precipitates (D.15 or C'.14 or D'10) with concentrated HNO<sub>3</sub>. Add 5 mL of H<sub>2</sub>O<sub>2</sub> to facilitate dissolution (possibly heat the sample if poor dissolution). Pour the solution into the annotated beaker. Rinse the centrifuge tubes with 8M HNO<sub>3</sub>.*

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*Adjust the pH to 1.5 with concentrated NH<sub>4</sub>OH. Close to pH = 1, a white precipitate appears. If the precipitate is not visible, add 3.5 mL of the CaCl<sub>2</sub> solution (100 mg.mL<sup>-1</sup>).*

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*Stir the solution for ½ hour while heating at 70 - 80 °C. Extend the heating by ½ hour.*

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*Rinse the centrifugation pellet with 50 mL of an aqueous solution of oxalic acid adjusted to pH 1.5.*

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*Redissolve the pellet with 10 to 30 mL of concentrated HNO<sub>3</sub> and 5 to 15 mL of H<sub>2</sub>O<sub>2</sub>. Pour into a 100 mL flask and boil at reflux (with a watch glass) for 2 hours until the red fumes disappear.*

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*Remove the watch glass and bring to dryness*

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*Redissolve the dry residue in 40 ml of 8M HNO<sub>3</sub>.*

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***During column loads, all the solutions collected at the bottom of the column are retained until the result of counting the sample.***

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*Store the washing solutions. Then, replace the recovery flask at the bottom of the column with a 125 mL cylindro-conical flask.*

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*Weigh 579.6 mg of  $\text{NH}_4\text{I}$  per sample, dissolve in 1 mL of ultrapure water, add 39 mL of 12M HCl<sup>1</sup>.*

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*Elute Pu with 40 mL of the 0.1M HCl /  $\text{NH}_4\text{I}$  solution prepared as previously described.*

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*Rinse the walls of the separating funnel and the column with 5 mL of 12M HCl.*

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*Add 5 mL of concentrated  $\text{HNO}_3$  to the column eluate and place the cylindro-conical flask on the sand bath. Bring to dryness at a gentle temperature (60 °C).*

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*Continue the treatment of the Pu fraction, by the steps G, L and M and the treatment of the Am-Cm fraction by the steps H, I, J, K, L and M.*

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*For a group of 4 samples, weigh 2.90 g of  $\text{NH}_4\text{I}$ . Dissolve in about 2 mL of ultrapure water. Transfer to a 10 mL graduated cylinder and make up to 5 mL. Transfer the solution to a 5 mL beaker. To elute 40 mL, add 1 mL of this solution to a 50 mL beaker and add 39 mL of 12M HCl. Pour into the separating funnel.*

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*Dissolve the residue from the cylindro-conical flask (step F.17) with concentrated HNO<sub>3</sub> and a few drops of H<sub>2</sub>O<sub>2</sub>.*

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*Bring to dryness and take up the dry residue with 10 mL of 8M HNO<sub>3</sub>*

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*Add 0.1 mL of 130 V H<sub>2</sub>O<sub>2</sub>. Cover with a watch glass and heat until the end of the effervescence to flush out excess H<sub>2</sub>O<sub>2</sub>.*

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*If purple vapors appear or if a pink-violet color of the solution persists, there is still iodine: add H<sub>2</sub>O<sub>2</sub> by fraction of 0.1 mL by heating at reflux until complete disappearance of the color.*

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*Add about 100 to 200 mg of NaNO<sub>2</sub>. Let stand at least ½ h.*

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*Prepare a column of 8 mm inside diameter with **1.5 g** of **AG1X8** resin (see Appendix A1 Preparation of columns).*

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*Put a recovery flask at the end of the column.*

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*Condition the column with 25 mL HNO<sub>3</sub> 8 M. (Discard the collected solution)*

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*Place the small flow tip with tap. Slow the flow at the bulb, 1 drop per sec.*

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*Pour the sample into the separating funnel and pass it on the column.*

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*Then pass 3 x 10 mL of 8M HNO<sub>3</sub> on the column by flushing the cylindro-conical flask that initially contained the sample.*

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***Continue processing of the Pu fraction by step L, source preparation step***

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*For a group of 4 samples, weigh 1.16 g of NH<sub>4</sub>I for a volume of 5 mL (For a group of 6 samples weigh 1.62 g for a volume of 7 mL). Dissolve them in about 2 mL of ultra-pure water. Transfer to a 10 mL graduated cylinder and make up to volume (5 mL for a group of 4, 7 mL for a group of 7). Transfer the solution to a 10 mL beaker. Introduce 1 mL of this solution into a 25 mL beaker and add 15 mL of 12M HCl.*

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*Take up the rust residue (Am-Cm fraction, Erlen F.9) with 10 mL of 9M HCl.*

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*Add 2 drops of H<sub>2</sub>O<sub>2</sub> 130 V. Cover with a watch glass.*

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*Heat for a few minutes (gentle temperature 60 °C). Cool without the watch glass.*

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*Prepare a double column with an inside diameter of 10 mm with **6.5 g of AG1X8** resin and **2 g of AG50W** resin (see Annex A1 Preparation of columns).*

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*Install a recovery flask at the outlet of the column.*

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*Condition the column with 30 mL of 9M HCl. Discard the bottom of the column.*

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*Place a 250 mL beaker under the column. Place a high flow tip at the bottom of the column.*

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*Load the solution into the separatory funnel and allow it to percolate on the column (Am-Cm not fixed).*

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*Rinse with 2 × 10 mL then 1 × 15 mL of 9M HCl by rinsing the Erlenmeyer flask initially containing the sample and collect into the beaker.*

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### *Precipitation of aluminium hydroxide*

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*Take the eluate at the column outlet (step H.9)*

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*Add 5 mL of 1M  $\text{Al}(\text{NO}_3)_3$ . Adjust the volume of the solution to about 150 mL with ultra-pure water.*

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*Adjust the pH to 8 with concentrated  $\text{NH}_4\text{OH}$ .*

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*Dissolve the pellet with 10 mL of concentrated  $\text{HNO}_3$  and 2 mL of  $\text{H}_2\text{O}_2$ . Reflux for 15 minutes. Pour into a 100 mL beaker and dry to 70 - 80 °C*

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*For samples containing few metals (herbs, small sample sizes, ...), continue the treatment with steps J.1 to J.7 and then step K. For samples rich in metals (sediments, soils, etc.), continue the treatment with steps J.8 to J.14 and then step K*

**J. Séparation Am / métaux Cu, Ni : Colonne TRU en milieu  $\text{HNO}_3$  2M,  $\text{Al}(\text{NO}_3)_3$  0,5M**

***Separation Am / métaux Cu, Ni : TRU Column in 2M  $\text{HNO}_3$ , 0,5M  $\text{Al}(\text{NO}_3)_3$  medium***

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*Dissolve the residue with 2.5 mL of 8M  $\text{HNO}_3$ . Heat to reflux until dissolved. Remove from the heat source and add a tip of spatula of ascorbic acid.*

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*Add 7.5 mL of ultrapure water to obtain 10 mL of 2M  $\text{HNO}_3$ - 0.5M  $\text{Al}(\text{NO}_3)_3$ .*

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*If the residue persists, filter the sample on a cellulose acetate filter (or pre-filter system installed on a syringe) before passing it through the column.*

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*Prepare a disposable 8 mm ID column with 0.70 g of TRU resin. For proper handling, it is recommended to wet the resin directly with 2M  $\text{HNO}_3$  before pouring it into the column (see A1 Preparing the columns) and putting glass wool on top of the resin.*

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*Condition the column with 10 mL of 2M  $\text{HNO}_3$ .*

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*Place a recovery flask to keep until the result. Load the sample on the column by volume of 2 ml, then rinse with  $2 \times 5$  mL of 2M  $\text{HNO}_3$ .*

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*Collect the eluate in a small 25 mL erlen.*

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*Dissolve the residue with 7.5 mL of 8M  $\text{HNO}_3$ . Heat slightly under reflux to dissolve the residue. Remove from the heat source and add a spatula tip of ascorbic acid.*

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*Add 10 mL of 1M  $\text{Al}(\text{NO}_3)_3$  and 12.5 mL of ultrapure water to dissolve the sample in 30 mL of 2M  $\text{HNO}_3$  and 0.5M  $\text{Al}(\text{NO}_3)_3$ .*

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*If the residue persists, filter the sample on a cellulose acetate filter (or pre-filter system installed on a syringe) before loading it through the column.*

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*Prepare a 8 mm ID column with 2.0 g of TRU resin. For easier handling, it is recommended to wet the resin directly with 2M  $\text{HNO}_3$  before pouring it into the column (see A1 Preparing the columns) and putting glass wool on top of the resin*

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*Condition the column with 30 mL of 2M  $\text{HNO}_3$ .*

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*Place a recovery flask to keep until the result. Load the sample on the column, then rinse with  $2 \times 5$  mL of 2M  $\text{HNO}_3$ .*

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*Elute with 40 mL of 0.03M  $\text{HNO}_3$ . Collect the eluate in a small 50 mL erlen.*

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*Reduce the eluated solution of the TRU column (step J.7 or J.14) to dryness at a gentle temperature (60 °C).*

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*Dissolve the dry residue to obtain a mixture of 20 mL of 1M  $\text{HNO}_3$  / 93%  $\text{CH}_3\text{OH}$ .*

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*Prepare a column of 8 mm inner diameter (ID) with 4 g of AG1-X4 resin (see Annex A1 Preparation of columns). For easier handling, wet the resin directly with 10 mL 1M  $\text{HNO}_3$  / 93%  $\text{CH}_3\text{OH}$  before transferring to the column.*

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de  $\text{HNO}_3$  1M/ $\text{CH}_3\text{OH}$  93 %.

*Remove the air bubbles trapped in the resin with a very fine glass rod after having previously closed the tap, adding if necessary up to 25 mL of 1M  $\text{HNO}_3$  / 93%  $\text{CH}_3\text{OH}$ .*

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*Before the end of the conditioning, remove the small bubbles using a very fine glass rod after having previously closed the tap. **The column should never be dry.***

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*Load the sample on the column. Monitor the liquid level. **The column should never be dry.***

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*Elute the rare earths with 80 mL of 0.1M  $\text{HCl}$  / 0.5M  $\text{NH}_4\text{SCN}$  / 80%  $\text{CH}_3\text{OH}$ , prepared at the last minute.*

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*Place a 100 mL Erlenmeyer at the bottom of the column.*

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*Bring the eluated solution to dryness. Dissolve the residue with concentrated  $\text{HNO}_3$  +  $\text{H}_2\text{O}_2$  and bring back to dryness (at a gentle temperature of 60 °C).*

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*Add a few mL of concentrated  $\text{HNO}_3$  then 2.5 mL of 5%  $\text{HNaSO}_4$ , and bring back to dryness*

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*Continue processing the Am-Cm fraction by step L, source preparation step*

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*Dissolve the dry residue (step G.17, Pu and K.15, Am) in 0.2 mL of concentrated HNO<sub>3</sub> and 2.5 mL of ultra-pure water. Transfer to an annotated glass liquid scintillation vial (Pu or Am fraction). Rinse the cylindro-conical flask or Erlenmeyer flask with twice 3.5 mL of 15% Na<sub>2</sub>SO<sub>4</sub> and then 0.5 mL of deionized water (add the rinse solution to the scintillation vial).*

- L.2** Ajouter 5 gouttes de NH<sub>4</sub>OH concentré. Add 5 drops of concentrated NH<sub>4</sub>OH.

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*Adjust the pH of the solution to 1.8 - 2.0 with the pH meter equipped with a previously calibrated microelectrode. Rinse the electrode with 0.5 mL of ultra-pure water above the scintillation vial.*

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*If it is not a mirror-polished disc, etch the stainless steel disc for 15 minutes at 1.2 A with 10 mL of the pickling solution. Discard the pickling solution. Rinse the system repeatedly with ultra-pure water.*

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*Place the disc in the disc holder and screw on the plastic bottle. Check tightness with ultra-pure water.*

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*Transfer the solution to be electrolyzed into the plastic bottle. Rinse the liquid scintillation vial with 0.5 mL of ultrapure water.*

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*Lower the platinum wire (anode) through a funnel (the end of which has been cut) until the wire is in contact with the stainless steel disc and then raise the anode by 5 mm.*

*Electroplate the solution for 2 hours at 1 A, with an anode rotation speed of about 40 rpm. Stop the electrolysis by adding 1 ml of concentrated NH<sub>4</sub>OH.*

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*Raise the platinum wire out of the vial and then only turn off the power.*

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*Remove the stainless steel disc and the plastic bottle containing the rest of the solution. Pour the solution into the liquid scintillation vial initially containing the solution to be electrodeposited and store until the results are validate.*

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*Extract the stainless steel disc. Rinse with ultra pure water and ethanol.*

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*Resorb a large portion of the liquid by tilting the edge of the stainless steel disc on absorbent paper. Place the stainless steel disk on the stainless steel drying rack. Dry the discs by placing the assembly on a hot plate (250 to 300 °C).*

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*On the back of the disc, identify the sample/fraction.*

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*If an alpha chamber is available, start counting\* immediately. Otherwise, keep the stainless steel disc in a plastic box with a lid until it is counted.*

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*Remove the platinum wire from the equipment and install it in the corresponding beaker containing 8M HNO<sub>3</sub>. Let it heat ¼ hour.*

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*Rinse this platinum wire with ultra-pure water.*

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\*Counting. In LMRE Pu and Am-Cm isotopes are determined with a low level background alpha spectrometer (alpha-analyst, Canberra with Genie 2000 software). Detectors are 450 mm<sup>2</sup> PIPS detectors, with an efficiency around 34% for <sup>239</sup>Pu. Counting time was 14 days

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#### **Preparation of a column made of a single resin**

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#### **Preparation of a column made of a two resins**

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## Annex 2. PRESENCE D'INTERFERENCES SPECTRALES

### PRESENCE OF SPECTRAL INTERFERENCES

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*In the case of a spectrum showing interferences, carry out a new purification:*

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##### ***Dissolve the electroplating for that :***

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*Take off the label behind the source (if any).*

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*Remove from the heat source and add 3-4 drops of concentrated HCl.*

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*Rinse the stainless steel disk with this solution until complete recovery of the deposit (use a disposable pipette and disposable tweezers).*

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
*Restart the procedure at step G (Pu) or K (Am).*


## **4 REFERENCES**

Partly : F. Goutelard, M. Morello, D. Calmet (1998), Alpha-spectrometry measurements of americium and curium at trace levels in environmental samples using extraction chromatography. *Journal of Alloys and Compounds*, 271-273, pp.25-30.

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**F1. Diagram of the procedure for the measurement of the isotopes of Pu, Am and Cm**