RiO5 METHOD (12-13)

IRSN/LMRE/MESURE/H-3/P4010.12V15

IRSN/LMRE/MESURE/H-3/MO4010.12.1V9

Sabine Charmasson

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**Tritium (Free and Organically bound tritium) —Liquid Scintillation
— environmental samples**

Disclaimer

It is the responsibility of each analyst to follow established practices when handling and examining the samples referenced in this Rio5 Cookbook. Although the methods may have been tested by each laboratory identified as the source, each user must perform a validation procedure to ensure the validity of their results. Woods Hole Oceanographic Institution, its officers, directors and employees are not responsible for any of the data or the results that may be achieved from using the information in the Rio5 Cookbook and disclaim all liability for the same.

Table of Contents

1	SCOPE	1
2	EQUIPMENT AND CHEMICAL REAGENTS	1
2.1	EQUIPMENT	1
2.2	TRACERS/RADIOACTIVE SOURCES	2
2.3	CHEMICAL REAGENTS	2
2.4	SOLUTIONS	2
3	PROCEDURE	2
3.1	COMBUSTION POUR MESURE DU TRITIUM LIE A LA MATIERE ORGANIQUE (TOL)	3
3.1.1	PRINCIPE/ <i>PRINCIPLE</i>	3
3.1.2	PRÉPARATION DU MATERIEL / PREPARATION OF THE MATERIAL	5
3.1.3	COMBUSTION	7
3.2	MESURE/MEASUREMENT	9
4	REFERENCES	10
5	FLOW CHART	11
6	IMAGES	13

1 SCOPE

This method both in French (original language) and English (*italic characters*) specifies the minimum requirements and laboratory methods for the preparation of samples and material to carry out analysis of Tritium in environmental samples with liquid scintillation.

After collection in the field, the samples are either frozen or keep fresh until return to the laboratory taking the necessary precautions to avoid any exchange with the air (samples must be put in sealed bags under vacuum as soon as possible after sampling).

In order to determine tritium concentrations, environmental samples are generally first freeze-dried. The tritium obtained in the condensate is the HTO fraction of the sample (tissue free water tritium). Then the residue is burnt in a special tubular furnace with oxygen in the presence of a catalyst. The water fraction trapped, designated as Organically Bound Tritium (OBT), is adjusted at pH 7, distilled and, measured using a liquid scintillation counter (Pointurier et al., 2004; Cossonnet et al., 2009).

As the combustion phase is critical for OBT measurement, this document details the procedure for this step.

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The measurement of ^3H is carried out using liquid scintillation counters (PACKARD, QUANTULUS).

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See text of Procedure for other materials

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Most of the equipment used requires a particular maintenance characteristic of its operation. The result of these operations is reported in the device's life record.

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In the absence of any indication on their certificate, the lifetime of the sources is fixed at most half of the tritium period, i.e. 6 years

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3 PROCEDURE

Warning

The original documents quote various internal documents relating in particular to the calibration, the maintenance of the equipment used as well as aspects of internal quality assurance in the laboratory. These documents are not published here and it is appropriate for each laboratory that uses these procedures to verify that it works in the state of the art.

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Tritium is extracted from the samples as tritiated water and then measured by liquid scintillation. Water samples (freshwater, rain, drinking water) are measured directly or after distillation (generally seawater measurement is not possible by 'direct' liquid scintillation due to very low amount of tritium in the marine compartment. It thus requires the use of other methods such as liquid scintillation after electrolytic enrichment or gas mass spectrometry by measuring the ^3He grown in from tritium decay). The other matrices undergo a transformation (freeze-drying, combustion, distillation) which makes it possible to extract the hydrogen which they contain in the H_2O form.

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The procedure is shown schematically in FIG. F1, depending on the nature of the sample.

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The quantities of water collected at each stage are measured by weighing. The mass / volume correspondence is performed using a density of 1 g.mL^{-1}

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For organically bound tritium, the analysis of each sample may be repeated once or twice.

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Procedure of the combustion using a tubular furnace (Pyroxydizer) of the samples, in order to recover hydrogen bound to organic matter in the form of water (MO4010.12.1).

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
Each combustion line consists of two heating units (furnace F1 and furnace F2)

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
In each line is arranged a quartz tube. In the part of the tube located at the furnace F1 is the combustion boat containing the sample and in that at the furnace F2 is the catalyst (oxidizing mixture alumina / Pt / crushed quartz) (Figure F2).

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The furnace F1 is semi-insulated and ventilated, the furnace F2 is totally insulated

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The purpose of the combustion is to recover water (containing tritium) by catalytic oxidation of the combustion gases. The combustion takes place in two stages:

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An argon stream sweeps the sample while a small stream of oxygen scavenges the catalyst. The argon flow drives the combustion gases. The flow of oxygen is kept low enough to prevent combustion from heating up

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an oxygen flow sweeps the sample and the catalyst for total oxidation of the product gases, while a small argon flow sweeps the sample.

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The produced water is then cold trapped (Figure F2).

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The operation of the furnaces is controlled by the software of the supplier. The use of thermocouples allows temperature control.

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The catalyst, an alumina / platinum mixture, is introduced into the tube by means of a funnel. Quartz wool (or glass wool) is placed on either side of the catalyst to keep it locked.

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Before any combustion, the mixture is heated.

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After one or more combustion of samples, the quality of the catalyst is checked visually; he must be gray. If it is white, it must be renewed.

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The heating of the tube without sample must be realized :

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either after a prolonged stop of the installation, to eliminate any trace of humidity,

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or after finding the presence of black residues (poor combustion and as mentioned above after any renewal of catalyst).

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

After weighing the 50 ml round flask and the U-tube, assembly is mounted. The round flask and the U-tube dive into a Dewar filled with dry ice. It is advisable to use high form Dewar to trap the water properly. The bottom of the U-tube should just be in contact with the dry ice.

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F3 : Assembly for the collection of water after combustion

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The empty combustion boat is weighed.

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A given mass of sample is added. This mass is a function of the available sample mass, the nature of the sample and the percentage of hydrogen to collect, if possible, at least 10 g of water (see guidance § 3.1.3.5).

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The combustion boat is introduced into the quartz tube, at the furnace F1, in abutment with the quartz wool (Figure F2).

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The time between the weighing of the sample and the introduction of the combustion boat into the furnace should be as short as possible to avoid any rehydration of the sample.

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The programming is designed so that the preheating of the furnace F2 at 800 °C is started, the temperature of the furnace F1 being maintained at about 35 °C and the furnace ventilated.

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This step carried out, the combustion sequence of the Furnace F1 starts up to reach 850°C and the temperature of the furnace F2 is raised to 1000 °C.

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The cryostat can be programmed according to the combustion sequences and it is also controlled with the software.

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It is advisable to heat the water of the cryostats 2 hours before the end of the sequence to dry the glassware above the round flask and the U-tube.

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When there are no more traces of water in the assembly :

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When, at the end of the combustion, the water collected is slightly colored (yellow/orange), a 2nd combustion of the totality of the water is carried out to purify it.

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Before the measurement by scintillation, the water collected is neutralized to pH 7 if the pH is acidic with addition of dry CaCO₃, then distilled.

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Combustion blanks are regularly realized each year in order to verify that there is no tritium contamination in the equipment (About 12 g of Abatilles water (deep spring water) are burned and the collected water is analysed).

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Indicative quantity and approximate duration of combustion according to the nature of the samples.

Echantillon/Sample	Quantité (g) Quantity (g)	Durée totale Total duration
Crustacés/ <i>Crustaceans</i> Mollusques/ <i>Molluscs</i> Mousses aquatiques/ <i>Aquatic mosses</i> Phanérogames/ <i>marine phanerogams</i>	30 – 40	6h00
Matrices sucrées (miel)/ <i>Sweet matrices (honey)</i> Poissons gras (anguille)/ <i>Fatty fish (eel)</i> Fromages gras/ <i>fatty cheese</i>	15 – 20	6h00
Herbe, mousse terrestre, feuille d'arbre, bois <i>Grass, moss, tree leaf, wood</i> Fruit, légume, salade, champignon, céréale <i>Fruit, vegetable, salad, mushroom, cereal</i> Poisson non gras (chair) <i>Non-fatty fish (flesh)</i>	30 - 40	6h00
Eau (eau colorée après 1 ^{ère} combustion) <i>Water (colored water after 1st combustion)</i>	Eau recueillie	4h00
Blanc de combustion/ <i>Combustion blank</i>	≥ 12	2h30

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The tritium activity in the water is determined by measuring the β^- emission of ^3H by liquid scintillation

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Samples can be associated in a group for counting. In IRSN-LMRE to each counting group are inserted:

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According to the manufacturers' recommendations, the counters are checked using 2 vials of sealed sources (^3H , ^{14}C), bearing a label indicating their date of manufacture and the activity they contain and a sealed vial for the background noise (Background label, date of manufacture).

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These checks are performed for maintenance purposes (for example, threshold values given by the manufacturer for counter efficiencies for ^3H , 56% for Packard and 55% for Quantulus) and are used by the person (supplier) who performs maintenance or intervention. These data are archived on the computer associated with the device.

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Activities are always calculated on the date requested by the customer in the test request or, failing that, on the date of sampling, taking into account the tritium period. The results are always expressed in Bq.L^{-1} . In the case of organically bound tritium, the values in Bq.kg^{-1} sec are given for information only.

4 REFERENCES

Pointurier, F., Baglan, N., Alanic, G., Chiappini, R., 2004. A method for the determination of low-level organic-bound tritium activities in environmental samples. Applied Radiation and Isotopes 61, 293–298.

C. Cossonnet C., NeivaMarques A.M., Gurriaran R., 2009: Experience acquired on environmental sample combustion for organically bound tritium measurement. Applied Radiation and Isotopes 67, 809–811

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NF ISO 9698 (February 2011) « Water quality – Determination of tritium activity concentration – liquid scintillation counting method »

XP M 60-824 (AFNOR, July 2016) TFWT and OBT analytical procedure in environmental samples

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
Figure F1: procedure for the analysis of ³H in the different samples (SL: liquid scintillation)

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