RiO5 METHOD (50)

LRC-AT-MOP-006-ind 1

Sabine Charmasson

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sabine.charmasson@irsn.fr

This is a laboratory operating protocol developed by the Radioecology Laboratory of Cherbourg, France at IRSN/PSE-ENV/SRTE

Contributor:

Pascal Bailly du Bois (pascal.bailly-du-bois@irsn.fr)

Cesium, Antimony, Ruthenium and Cobalt —Gamma spectroscopy— large volume water sample

Disclaimer

It is the responsibility of each analyst to follow established practices when handling and examining the samples referenced in this Rio5 Cookbook. Although the methods may have been tested by each laboratory identified as the source, each user must perform a validation procedure to ensure the validity of their results. Woods Hole Oceanographic Institution, its officers, directors and employees are not responsible for any of the data or the results that may be achieved from using the information in the Rio5 Cookbook and disclaim all liability for the same.

Table of Contents

1	SCOPE	1
2	EQUIPMENT CHEMICAL REAGENTS	1
2.1	EQUIPMENT	1
2.2	TRACERS	1
2.3	CHEMICAL REAGENTS	1
2.4	SOLUTIONS	2
3	PROCEDURE	2
4	REFERENCES	9
5	FLOW CHART	10
6	IMAGES	11

1 SCOPE

The method presented here specifies the minimum requirements and laboratory methods for preconcentration before gamma-spectrometry measurement of dissolved cesium, antimony, cobalt and ruthenium in large volume water sample.

It is based on the co-precipitation of MnO_2 + KCFC. MnO_2 co-precipitates the hydrolysed forms of cobalt, ruthenium and antimony. Cesium is co-precipitated by addition of fine particles of KCFC (<100 μm , 8 to 10 g). This precipitate, brought to pH 3 - 3.5, makes it possible to obtain the best extraction yields, in particular for antimony. In order to obtain good extraction yields, the precipitate is also kept in suspension for a minimum of 3 hours (stirring by air bubbling) before decanting (Gandon and Guéguéniat, 1992; Bailly-du-Bois et al., 1999)

2 EQUIPMENT AND CHEMICAL REAGENTS

2.1 Equipment

- for ex MILLIPORE YY30 293 16 natural water filtration system equipped with a pressure tank at 2.5 bar. The soluble phase (less than 0.45 μm) is isolated from suspended solids from 120 to 480 liters of sample;
- Containers for collecting concentrates of dissolved elements. 10 or 20 liters capacity
- Compressed air generator
- 120 liters reaction tank with bubble rod connected to the compressed air generator;
- field portable pH meter

CLEANING OF BARRELS

- The hydrated manganese dioxide and copper ferrocyanide attached to the tank/cubitainer and other material is facilitated when the reagents are in hydrated form. The manganese dioxide is solubilized as a Mn^{2+} ion by a reducing acid mixture (HNO_3 0.1N + H_2O_2 at 110 volumes diluted to 5%)

2.2 Tracers

- No tracer/ the extraction yield for antimony and cesium is about 97% and 100% respectively

2.3 Chemical reagents

- Potassium permanganate
- Hydrogen peroxide 110 volumes;
- KCFC (powder < 100 μm : ferrocyanide of cobalt and potassium)
- Nitric Acid 65%
- Sodium hydroxide

2.4 Solutions

- potassium Permanganate (solution 30g l⁻¹)
- Hydrogen peroxide 110 volumes;
- KCFC (ferrocyanide of cobalt and potassium) (10 g pre-weighed flask)
- Nitric Acid 65% and Sodium hydroxide 1N for pH adjustment

3 PROCEDURE

Prerequisite

The protocols relate to seawater samples with natural salt concentration up to 36 g.l⁻¹ and whose pH is 8.2. For fresh water samples with little dissolved salts, add sodium chloride after filtration to reach or exceed a salinity of 3-4 g.l⁻¹. in order to allow precipitation of the manganese dioxide and copper ferrocyanide.

It is also important to add sodium hydroxide to the filtered sample in order to be placed in a basic medium when the filtered water is at a pH below 7.


This method applies for a sample volume of 120 l but they can be applied to other volumes by adjusting the proportions of reagents.

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Fill a plastic container with 120 liters of seawater. As there is no risk of pollution, the in-board fire pump could be employed. In open seas with low quantities of suspended matter, it is not necessary to filter the seawater. Record at the same time the location and depth and put a label corresponding to the ID station number on the tank. Secure strongly the containers to allow processing for 6 hours.

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To obtain a good mixing of seawater during several hours, the easiest method is to introduce compressed air at the bottom of the tank with thin bubble rods. Several tanks could be mixed at the same time connecting them to the compressed air generator with a thick rubber vacuum tube where bubble rods and their flexible pipe are pricked by means of surgical needles. To stop mixing, just unplug the needles. Air pressure of 1-2 bar ($\text{kg}\cdot\text{cm}^{-2}$) is sufficient.

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*Chemical adsorbent is a solution of potassium permanganate (30 g KMnO_4 in 1 liter or 30g powder) for measurement of *Ru*, *Co*, *Sb*; mixed with a powder of double cobalt – potassium ferrocyanide ($\text{Fe}(\text{CN})_6\text{CoK}$, 10 g) for caesium measurement. This powder is insoluble, it must be fine enough to stay in suspension during agitation and ensure a good exchange with seawater.*

The water becomes deep purple after their introduction. If permanganate is introduced as powder, it must be slowly introduced in order to dissolve it properly. It is the same for ferrocyanide which must stay well mixed during hours and do not settle immediately.

3.4 Ajout H₂O₂

Add H₂O₂

Fig. 4 Transformation of potassium permanganate into manganese dioxide

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Within the tank KMnO₄ must be precipitated in manganese dioxide (MnO₂) to fix Sb, Ru, Co. The form obtained during the reduction of Mn (VII) in Mn (IV) by H₂O₂ in basic medium strongly depends on the operating protocol: quantity of H₂O₂ and speed with which this reagent is introduced, chronology in the addition of reagents.

In order to obtain more reproducible forms:

- if the KMnO₄ is added as a solution, pour 100 mL of H₂O₂ **BEFORE** adding the KMnO₄ solution,
- if the KMnO₄ is introduced as powder pour **AFTER** the powder **VERY SLOWLY** about 100 mL of concentrated hydrogen peroxide (H₂O₂ 110 vol).

The mixture goes into effervescence and becomes brown with abundant foam (the MnO₂ catalyzes the degradation of H₂O₂). If a minute later, the water remains slightly red-purple, add (by small additions) a little more H₂O₂ until the solution is true brown.

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*When the precipitate is formed, the pH of the medium has to be adjust to 3 - 3.5. Nevertheless, it is extremely important that all the H_2O_2 in excess has been eliminated **before**. Indeed, in an acid medium, H_2O_2 behaves as a reducing agent with respect to Mn (IV). The presence of H_2O_2 at this stage would solubilize a part of MnO_2 . In order to this to be achieved it is better to wait 10 min before adjusting the pH.*

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The adsorption ratio is more efficient at pH 3.5 (100% for antimony). The pH is lowered by addition of concentrated nitric acid (HNO_3) or ammonia (NH_4OH) after complete disappearance of bubbles.


Fig. 5 Adjust pH 3.5 with NITRIC ACID 65 %

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After these treatments, close the tank and maintain mixing with bubbles during at least 3 hours.



Fig. 6 3 Hours mixing (bubbling air)

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Remove the bubbling rod and allow the precipitate to settle (It could take 3 to 10 hours, depending of movements of the boat and salinity.)

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When decantation is finished, water is clear and there is a net difference between supernatant and precipitate. Siphon the supernatant from the top by using a flexible tube. At this time, it is recommended to save some water in order to rinse the tank at the end. When the water level is close to the precipitate, take care not to move it and stop before its aspiration.



Fig. 8 *Eliminate the supernatant by aspiration of water from the surface*

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Large filters are placed in funnels and identified with the station number; the precipitate is then collected inside. It is important to rinse properly the ferrocyanide powder at the bottom of the tank (it contains all the caesium). After complete recovery and draining, filters can be fold for stocked in plastic bags for transport.


Fig. 9 *Filtrate the precipitate with large filters and fold the filter*

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*All the steps described so far are carried out **in the field**. However in case of bad weather precipitates (previous step) can be kept in containers of 20 liters for later filtration in the laboratory.*

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The filter is dried in an oven (between 50 and 90 °C). The precipitate is then collected (see photos) in gamma counting geometry and ready for measurement.

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A solution of H₂O₂ and HNO₃ is used to clean the tanks between each sample (about 50 mL of each in 5 liters).

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Bailly duBois, P., Guéguéniat, P., 1999. Quantitative assessment of dissolved radiotracers in the English Channel: sources, average impact of la Hague reprocessing plant and conservative behaviour (1983, 1986, 1988 and 1994). Continental Shelf Research 19, 1977–2002.

Gandon, R., Guéguéniat, P., 1992. Preconcentration of ¹²⁵Sb onto MnO₂ from seawater samples for Gamma-ray spectrometric analysis. Radiochimica Acta 57, 159}164.

5 FLOW CHART

Use of protections to prevent injuries from hydrogen peroxide, nitric acid and sodium hydroxide (gloves, mask)

PRECONCENTRATION ON MANGANESE IV DIOXIDE OF SOLUBLE ELEMENTS IN NATURAL WATERS

120 LITRES filtered seawater (0.45 µm) Fig. 1

Mixing by bubbling air, the rod must reach the bottom of the tank to prevent sedimentation of the KCFC powder Fig. 2

(7 < pH < 8) (adjust pH if necessary)

ADD 8 to 10g of KCFC by dispersion on the water surface

Add 30g or 1l of POTASSIUM PERMANGANATE SOLUTION (30g.l⁻¹) Error! Reference source not found.

Mixing FOR 10 MINUTES (bubbling air)

Add slowly* hydrogen peroxide in order to obtain complete transformation of potassium permanganate into manganese dioxide (the purple- pink color must completely disappear, only brown color remain Fig. 4 Normally less than 100 ml is enough.

STIRRING FOR 10 MINUTES (bubbling air)

to eliminate H₂O₂ en excess

Adjust pH 3.5 with NITRIC ACID 65 % Fig. 5 Error! Reference source not found.

in order the chemical yield for antimony adsorption reach 100%

3 Hours STIRRING (bubbling air) Fig. 6

for ¹²⁵Sb et ¹⁰⁶Ru adsorption

5 Hours DECANTATION of MnO₂ Error! Reference source not found.

ELIMINATE THE SUPERNATANT by aspiration of water from the surface Error! Reference source not found.

FILTRATE THE PRECIPITATE with large filters Error! Reference source not found.

* It is advisable to first pour 25 ml, wait until the pH is constant, then pour GENTLY to pH 3-3.5. The water color should be brown if still violet add a LITTLE MORE hydrogen peroxide H₂O₂.

PRECONCENTRATION ON MANGANESE IV DIOXIDE OF SOLUBLE ELEMENTS IN NATURAL WATERS

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Recovery of the precipitate after drying